

VNS Therapy™ Parameters for Clinical Response in Epilepsy



Fahoum, F., Boffini, M., Kann, L., Faini, S., Gordon, C., Tzadok, M., & El Tahry, R. (2022). VNS parameters for clinical response in Epilepsy. Brain stimulation.

Key Take Away 1

Retrospective analysis of VNS Therapy™ parameters revealed a dose of VNS for epilepsy associated with seizure frequency reduction that lies within the target range in the product labeling. Patients titrated near this dose experienced greater seizure frequency reduction than those titrated above or below it.

	N	12-month Responder Rate	Median Seizure Reduction
Output Current <1 mA	44	36%	34.46%
Output Current 1.5-1.75 mA	392	47%	43.27%
Output Current >= 2.5 mA	32	41%	32.76%

Clinical outcomes of people with VNS titrated to settings near the model-selected dose of 1.625mA. Patients in each group were selected to have the listed output current at any pulse width at 12-months of follow up. The response rate was calculated at 12 months after implant.



A Generalized Linear Mixed Model (GLMM) was used to identify relationships between both VNS parameters and demographic features to clinical response (≥50% reduction in seizures). This is a population-level outcome, indicating that some patients will require VNS intensities above or below this value. The left figure demonstrates the population level mean and standard deviation of outcomes predicted by the GLMM. The right figure displays the model's predicted outcome for several hypothetical patients with specific demographic parameters.

Key Take Away 2

VNS responders titrated to the appropriate VNS parameters at or before 12 months after implant had a durable response to the therapy at future follow-up. Those who stayed at settings within the 1.5mA-1.75mA range had the highest likelihood of prolonged response to VNS.

	N (Unique Visits After 12 Months)	Responder Rate	Median Seizure Frequency Reduction
Output Current <1.5 mA	7	71%	100%
Output Current 1.5-1.75 mA	209	87%	86.13%
Output Current > 1.75 mA	98	80%	75.72%

Long-term clinical outcomes of VNS responders titrated to the target intensity of VNS (1.5mA – 1.75mA) by 12 months after implant (n = 184 subjects). Patients initially titrated to this dose may have increased or decreased their VNS dose at follow up visits after 12 months.

Study Summary

Objective

Significant variance in VNS Therapy™ parameter settings was detected in post-market surveillance of VNS Therapy management. The most recent VNS guidelines from the American Academy of Neurology recommended further investigation of VNS parameters could improve the quality of VNS management and patient outcomes¹. This analysis was undertaken with the purpose to provide an evidence-based justification of a target dose of VNS Therapy to improve epilepsy outcomes.

Method

A total of 1178 subjects from LivaNova-sponsored clinical studies were assessed for this work. A generalized linear mixed model was employed to assess the correlation between programming settings and demographic features and clinical response (defined as a 50% reduction in seizure frequency). This regression model was selected because it is known to be suitable when there is non-independence between observed values, such as when a single patient contributes multiple data points toward an outcome. The model included quadratic factors for output current and duty cycle, which were modeled as continuous variables, while pulse width and signal frequency were modeled as categorical variables due to a limited number of available settings.

Results

The mixed model successfully accounted for 86% of the variance in the database population, indicating a robust model fit with appropriate variable selection. A population level output current and duty cycle were identified as correlating with a peak level of response to VNS Therapy. The model did not find that any of the pulse widths or signal frequencies assessed had a significant impact on response.

While the outcome of this analysis differs slightly from the manufacturer's labeling, it can be said to generally confirm the target dose range recommended in the VNS Therapy System Epilepsy Physician's Manual, and patients who do not respond at the introductory duty cycle of 10% should have their duty cycle increased. This outcome was independent of a time effect, indicating that even patients who have been under-dosed for some time could still benefit from achieving the target dose. Post-hoc analysis further indicated that patients who were programmed to the target dose and kept there tended to stay in response (>85% likelihood).

Conclusion

This analysis supports wider adoption of current dosing recommendations, and specifically the use of target output currents within the range defined in the VNS Therapy System Epilepsy Physician's Manual. There are no robust data available at present to advocate for the use of frequencies other than 20, 25, or 30 Hz in epilepsy for the purpose of maximizing clinical response. While there appears to be a positive effect of duty cycle on seizure reduction, duty cycle should be increased only if patients do not respond at 10% duty in favor of battery savings.

Limitations

The principal limitation of this retrospective analysis is that it utilized data collected from a variety of clinical studies of VNS Therapy. The studies included interventional and observational designs, different follow-up durations, targeted patients of slightly different demographic profiles, and the methods for data collection were not uniform across all studies. None of these studies were prospectively designed for the purpose of assessing the relationship between VNS Therapy parameters and clinical response.



Click here or scan
the QR code to read
the full paper

¹Morris, G. L., Gloss, D., Buchhalter, J., Mack, K. J., Nickels, K., & Harden, C. (2013). Evidence-based guideline update: vagus nerve stimulation for the treatment of epilepsy: report of the Guideline Development Subcommittee of the American Academy of Neurology. *Neurology*, 81(16), 1453-1459.

