

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon
Tip, Plastic Tubing, Female Luer Lock Connector, 5.0 mm, CP-22005



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

ORDER CODE	CP-22005
SIZE	5.0 mm
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes