

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Basket
Tip, Stainless Steel Shaft, Female Luer Lock Connector, 12 Fr, CP-
12012



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

ORDER CODE	CP-12012
SIZE	12 Fr
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes