

# Venous Return Cannulae, Optiflow Venous, 29 Fr, V182-29



The Venous cannulae are intended to be used to cannulate the major venous vessels during cardiopulmonary surgery for periods of up to six hours. Optiflow Cannulae are provided with a one-piece swirled and grooved tip with multiple holes. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative

ORDER CODE	V182-29
SIZE	29 Fr
AVAILABLE COATED	Depending on country registration
EFFECTIVE LENGTH (CM)	39
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes