Venous Return Cannulae, Single Stage, 90 Tip with Suture Collar, Wirereinforced Tubing, 10 Fr, V132-10

The Pediatric cannulae are intended to be used in both the venous and the arterial line of the extracorportal

The Pediatric cannulae are intended to be used in both the venous and the arterial line of the extracorportal circuit during cardiopulmonary surgery for periods of up to six hours. Specifically, Venous Cannulae are designed to drain blood from the patient during cardio-pulmonary bypass procedures. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

ORDER CODE	V132-10
SIZE	10 Fr
AVAILABLE COATED	Depending on country registration
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes
CONNECTOR ACCEPTANCE	3/16""-1/4

Contact us Complete the form to get in touch with a LivaNova representative.



Contact Us