

# Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 7.0 mm, CP-22007



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

|                         |          |
|-------------------------|----------|
| <b>ORDER CODE</b>       | CP-22007 |
| <b>SIZE</b>             | 7.0 mm   |
| <b>AVAILABLE COATED</b> | no       |
| <b>QTY PER BOX</b>      | 10       |
| <b>LATEX-FREE</b>       | yes      |
| <b>NO DOP</b>           | yes      |