

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 8.0 mm, CP-22008



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

ORDER CODE	CP-22008
SIZE	8.0 mm
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes