

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Basket Tip, Stainless Steel Shaft, Female Luer Lock Connector, 10 Fr, CP-11010



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

<b>ORDER CODE</b>	CP-11010
<b>SIZE</b>	10 Fr
<b>AVAILABLE COATED</b>	no
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes