

# Venous Return Cannulae, Dual Stage, Wire-reinforced Tubing with Obturator and with Connector, 32/36 Fr, V900-284



The Venous cannulae are intended to be used to cannulate the major venous vessels during cardiopulmonary surgery for periods of up to six hours. Dual Stage Venous Return Cannulae are provided with a basket tip and proximal and distal wire reinforcement. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative

<b>ORDER CODE</b>	V900-284
<b>SIZE</b>	32/36 Fr
<b>AVAILABLE COATED</b>	Depending on country registration
<b>EFFECTIVE LENGTH (CM)</b>	47
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes