


Venous Return Cannulae, Single Stage, Bullet Tip with Suture Ring,
Wire-reinforced Tubing, 32 Fr, V900-02



The Venous cannulae are intended to be used to cannulate the major venous vessels during cardiopulmonary surgery for periods of up to six hours. Single Stage Venous Return Cannulae are provided with straight and right angled tip, with or without reinforced tubing. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative

ORDER CODE	V900-02
SIZE	32 Fr
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes
CONNECTOR ACCEPTANCE	3/8