

# Vents, Pre-bent Left Ventricle Catheters, PVC Bullet Tip, with Stylet, 19 Fr, W900-171



The vent catheters are intended to be used to vent blood of the left heart atrium or the left ventricle during cardiopulmonary surgery for periods of up to six hours. Depending on country registration, an optional phospholipid surface treatment may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative

|                              |                                   |
|------------------------------|-----------------------------------|
| <b>ORDER CODE</b>            | W900-171                          |
| <b>SIZE</b>                  | 19 Fr                             |
| <b>AVAILABLE COATED</b>      | Depending on country registration |
| <b>EFFECTIVE LENGTH (CM)</b> | 30                                |
| <b>QTY PER BOX</b>           | 10                                |
| <b>LATEX-FREE</b>            | yes                               |
| <b>NO DOP</b>                | yes                               |