


PATIENT INFORMATION CARD



Name of Patient: _____

Implantation Date: _____

In a medical emergency, call

Doctor: _____

Phone: _____

This person is implanted with a vagus nerve stimulator.
See reverse for emergency information.

Full safety information is available at VNSTherapy.com/manuals or by calling a number below.

To Stop Stimulation:

Hold or tape the VNS Therapy magnet over the left upper chest where the generator has been implanted. The Normal Mode cycle begins again once the magnet is removed.

In a Medical Emergency:

Contact a physician and call the local emergency number or go to the nearest emergency room. **LivaNova cannot provide medical care or medical advice.**

For Providers:

- ▶ [1.800.332.1375](tel:18003321375) (USA/Canada)
- ▶ [+1.281.228.7200](tel:+12812287200) (Worldwide)
- ▶ [+32.2.720.95.93](tel:+3227209593) (Europe/EMMEA/ASPAC)

IM-7601364-EPI

For quick information on
VNS Therapy and MRI safety
please scan the QR Code.

