PATIENT INFORMATION CARD

VNS Therapy®		Name of Patient: Implantation Date:	
Th	erapy"	Model/Serial Numbers	Generator:
	Name of the latest of the late	In a medical emergency, call Doctor: Phone:	
		This person is implanted with a vagus nerve stimulator. See reverse for emergency information.	

MR Conditional Device:

This person is implanted with a VNS Therapy System. MRI can be safely performed provided specific guidelines are followed. Patients should contact their physician before scheduling an MRI. The programming physician is required to program the device before and after an MRI. Providers: Full MRI safety information is available in the physician manual, available at VNSTherapy.com/manuals or by calling a number listed below.

To Stop Stimulation:

Hold or tape the VNS Therapy magnet over the left upper chest where the generator has been implanted. Stimulation will resume when the magnet is removed.

In a Medical Emergency:

Contact a physician and call the local emergency number or go to the nearest emergency room. **LivaNova cannot provide medical care or medical advice.**

For Providers:

- 1.800.332.1375 (USA/Canada)
- +1.281.228.7200 (Worldwide)
- +32.2.720.95.93 (Europe/EMMEA/ASPAC)