

Arterial Cannulae, Flexible, Curved Tip with Suture Ring, Non-reinforced  
Tubing, 21 Fr/ 7.0 mm, NA-5537



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-5537
SIZE	21 Fr/ 7.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes



**Contact us** Complete the form to get in touch with a LivaNova representative.

Contact Us