

Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring,  
Non-reinforced Tubing, 21 Fr/ 7.0 mm, NA-3427



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-3427
SIZE	21 Fr/ 7.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes