


Arterial Cannulae, Optiflow Arterial, Curved Tip, Wire- reinforced Tubing, 24 Fr/ 8.0 mm, A282-80C



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code A282-80C

Size 24 Fr/ 8.0 mm

Available Coated Depending on country registration

Effective Length (cm) 25

Qty Per Box 10

Latex-Free yes

No DOP yes