

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Conical Tip, Stainless Steel Shaft, 1/4"" Tubing Connector, 9 Fr, P608-30



The Coronary perfusion cannulae are intended to individually cannulate the coronary ostia and to selectively deliver cardioplegic solution during cardiopulmonary bypass surgery for periods up to six hours. Depending on country registration, an optional phospholipid surface treatment may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative

ORDER CODE	P608-30
SIZE	9 Fr
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us