

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Conical Tip, Stainless Steel Shaft , Female Luer Lock Connector, 15 Fr, P618-50



Livanova Antegrade Cardioplegia Coronary Ostia Perfusion Cannulae are intended for cardioplegia delivery into the coronary ostia. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	P618-50
SIZE	15 Fr
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes