


Cardioplegia Cannulae, Retrograde, Self-Inflating PVC Balloon, 14 Fr, RCS-11214



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

ORDER CODE	RCS-11214
SIZE	14 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	30.5
QTY PER BOX	5
LATEX-FREE	yes
CONNECTOR	Luer lock



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us