

# MICS Cannulae, Aortic Root Long Needle, 14 Ga (7 Fr), AR-17014



Livanova Aortic Root Long Needle is specifically designed for minimally invasive cardiac surgery in lateral thoracotomy, especially mitral valve procedures. It's longer than a conventional aortic root (35 cm) which allows it to be positioned in lateral thoracotomy.

ORDER CODE	AR-17014
SIZE	14 Ga (7 Fr)
INSERTION NEEDLE	16 Ga- White
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	31,8
QTY PER BOX	10
LATEX-FREE	yes