

MICS Cannulae, Easy Flow DUO Femoral Arterial, 19 Fr, 103-210



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

Order Code 103-210

Size 19 Fr

Available Coated no

Effective Length (cm) 28.5

Qty Per Box 5

Latex-Free yes