

# Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 20 Fr, SU-12502



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

Order Code SU-12502

Size 20 Fr

Available Coated no

Effective Length (cm) 36,8

Qty Per Box 10

Latex-Free yes