

Thoracic Catheters, Right Angle, 12 Fr, TC-11112



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2''' increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-11112
SIZE	12 Fr
EFFECTIVE LENGTH (CM)	51,8
LATEX-FREE	yes
CONNECTOR ACCEPTANCE	3/16 in



Contact us Complete the form to get in touch with a LivaNova representative.

Contact Us