

Venous Return Cannulae, Dual Stage, Wire-reinforced Tubing, 36/50 Fr, RDS-61150



The Dual Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The Dual Stage Cannula has a proximal flexible polyvinyl chloride tubing segment for drainage of the right atrium. The catheter body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava. The proximal and distal segments of the catheter are available in wire or non-wire reinforced tubing.

Order Code	RDS-61150
Size	36/50 Fr
Available Coated	no
Effective Length (cm)	41
Qty Per Box	10
Latex-Free	yes
No DOP	yes