

Venous Return Cannulae, Single Stage, Straight Lighthouse Tip, Wire- reinforced Tubing, 26 Fr, RV-40026



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

Order Code RV-40026

Qty Per Box 10

Latex-Free yes