

Vents, Aortic Root Aspirating Needle, 16 Ga, AX-30321



The Air Aspirator is indicated for use in delivery of cardioplegia solution during cardiopulmonary bypass surgery and aspirate air from the aorta at the end of the cardiopulmonary bypass procedure

Order Code	AX-30321
------------	----------

Size	16 Ga
------	-------

Available Coated	no
------------------	----

Effective Length (cm)	5,1
-----------------------	-----

Qty Per Box	20
-------------	----

Latex-Free	yes
------------	-----