

Accessories, AB1534



Livanova Accessories include porous vent plugs, vented connector caps, rapid prime sets, tubing organizers, vasco loops, tourniquets and pressure monitor lines.

| | |
|------------|--------|
| Order Code | AB1534 |
|------------|--------|

| | |
|-----------------------|----|
| Effective Length (cm) | 40 |
|-----------------------|----|

| | |
|-------------|----|
| Qty Per Box | 50 |
|-------------|----|

| | |
|------------|-----|
| Latex-Free | yes |
|------------|-----|