

## Accessories, Pressure Monitor Lines, AB1517



Livanova Accessories include porous vent plugs, vented connector caps, rapid prime sets, tubing organizers, vasco loops, tourniquets and pressure monitor lines.

|                              |        |
|------------------------------|--------|
| <b>ORDER CODE</b>            | AB1517 |
| <b>EFFECTIVE LENGTH (CM)</b> | 215 cm |
| <b>QTY PER BOX</b>           | 20     |
| <b>LATEX-FREE</b>            | yes    |