

# Accessories, Pressure Monitor Lines, AB1517



Livanova Accessories include porous vent plugs, vented connector caps, rapid prime sets, tubing organizers, vasco loops, tourniquets and pressure monitor lines.

|            |        |
|------------|--------|
| Order Code | AB1517 |
|------------|--------|

|                       |        |
|-----------------------|--------|
| Effective Length (cm) | 215 cm |
|-----------------------|--------|

|             |    |
|-------------|----|
| Qty Per Box | 20 |
|-------------|----|

|            |     |
|------------|-----|
| Latex-Free | yes |
|------------|-----|