

# Arterial Cannulae, Flexible, Curved Tip with Suture Ring, Non-reinforced Tubing, 21 Fr/ 7.0 mm, NA-5517



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

<b>ORDER CODE</b>	NA-5517
<b>SIZE</b>	21 Fr/ 7.0 mm
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	25
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes