

Arterial Cannulae, Flexible, Curved Tip with Suture Ring, Non-reinforced
Tubing, 24 Fr/ 8.0 mm, NA-5518



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-5518
SIZE	24 Fr/ 8.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes