

# Arterial Cannulae, Flexible, Curved Tip with Suture Ring, Non-reinforced Tubing, 18 Fr/ 6.0 mm, NA-5526



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

<b>ORDER CODE</b>	NA-5526
<b>SIZE</b>	18 Fr/ 6.0 mm
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	25
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes