

Arterial Cannulae, Flexible, Curved Tip with Suture Ring, Non-reinforced
Tubing, 24 Fr/ 8.0 mm, NA-5538



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

| | |
|-----------------------|---------------|
| ORDER CODE | NA-5538 |
| SIZE | 24 Fr/ 8.0 mm |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 25 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |