

Arterial Cannulae, Flexible, Straight Tip with Moveable Suture Ring,
Non-reinforced Tubing, 18 Fr/ 6.0 mm, NA-3426



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

ORDER CODE	NA-3426
SIZE	18 Fr/ 6.0 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	25
QTY PER BOX	10
LATEX-FREE	yes