

# Arterial Cannulae, Flexible, Straight Tip with Suture Ring, Non-reinforced Tubing, 24 Fr/ 8.0 mm, NA-4528



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

<b>ORDER CODE</b>	NA-4528
<b>SIZE</b>	24 Fr/ 8.0 mm
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	25
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes