

# Arterial Cannulae, Flexible, Straight Tip with Suture Ring, Non-reinforced Tubing, 18 Fr/ 6.0 mm, NA-4536



The Aortic Arch Cannula is indicated for use in perfusion of the ascending aorta during cardiopulmonary bypass surgery.

|                              |               |
|------------------------------|---------------|
| <b>ORDER CODE</b>            | NA-4536       |
| <b>SIZE</b>                  | 18 Fr/ 6.0 mm |
| <b>AVAILABLE COATED</b>      | no            |
| <b>EFFECTIVE LENGTH (CM)</b> | 25            |
| <b>QTY PER BOX</b>           | 10            |
| <b>LATEX-FREE</b>            | yes           |