

Arterial Cannulae, PureFlex, Curved Tip, With Suture Flange, Wire-reinforced Tubing, 22 Fr/ 7.3 mm, RA22-2BC



The Arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary surgery for periods of up to six hours. The Pureflex Arterial cannulae consist of flexible tubing (tapered to a curved or straight flow tip) made of PVC with wire reinforcement

ORDER CODE	RA22-2BC
SIZE	22 Fr/ 7.3 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	21.6
QTY PER BOX	10
LATEX-FREE	yes