

# Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 4.0 mm, CP-22004



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

Order Code      CP-22004

Size                4.0 mm

Available Coated   no

Qty Per Box        10

Latex-Free          yes

No DOP            yes