

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Basket  
Tip, Stainless Steel Shaft, Female Luer Lock Connector, 10 Fr, CP-  
11010



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

ORDER CODE	CP-11010
SIZE	10 Fr
AVAILABLE COATED	no
QTY PER BOX	10
LATEX-FREE	yes
NO DOP	yes