

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Basket Tip, Stainless Steel Shaft, Female Luer Lock Connector, 14 Fr, CP-11014



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

|                         |          |
|-------------------------|----------|
| <b>ORDER CODE</b>       | CP-11014 |
| <b>SIZE</b>             | 14 Fr    |
| <b>AVAILABLE COATED</b> | no       |
| <b>QTY PER BOX</b>      | 10       |
| <b>LATEX-FREE</b>       | yes      |
| <b>NO DOP</b>           | yes      |