

## MICS Cannulae, Bi-Flow Femoral Arterial, 19 Fr, 203-019



Bi-Flow is the only bidirectional arterial cannula on the market designed to reduce the risk of limb ischemia in an easy, safe and reproducible way.

Unlike traditional cannulae, Bi-Flow's patented design incorporates a unique elbow and downstream perfusion channel. The cannula is easy to insert and remove and provides simultaneous systemic and distal perfusion to the lower limb.

<b>ORDER CODE</b>	203-019
<b>SIZE</b>	19 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	26,7
<b>QTY PER BOX</b>	5
<b>LATEX-FREE</b>	yes