

# MICS Cannulae, Bi-Flow Femoral Arterial, 19 Fr, 203-019



Bi-Flow is the only bidirectional arterial cannula on the market designed to reduce the risk of limb ischemia in an easy, safe and reproducible way. Unlike traditional cannulae, Bi-Flow's patented design incorporates a unique elbow and downstream perfusion

channel. The cannula is easy to insert and remove and provides simultaneous systemic and distal perfusion to the lower limb.

Order Code	203-019
Size	19 Fr
Available Coated	no
Effective Length (cm)	26,7
Qty Per Box	5
Latex-Free	yes