

## MICS Cannulae, Easy Flow DUO Femoral Arterial, 19 Fr, 103-210



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

<b>ORDER CODE</b>	103-210
<b>SIZE</b>	19 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	28,5
<b>QTY PER BOX</b>	5
<b>LATEX-FREE</b>	yes