

MICS Cannulae, Easy Flow DUO Femoral Arterial, 23 Fr, 103-310



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

| | |
|-----------------------|---------|
| Order Code | 103-310 |
| Size | 23 Fr |
| Available Coated | no |
| Effective Length (cm) | 28.5 |
| Qty Per Box | 5 |
| Latex-Free | yes |