

# MICS Cannulae, Stand-Alone Guidewire, 0.035 (in), IQ35F180J6F



Livanova Vascular Dilator Kit, for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035 x 180 cm guide-wire, a blade scalpel and 8,12,16,20 and 24 Fr dilators.

<b>ORDER CODE</b>	IQ35F180J6F
<b>SIZE</b>	0.035 (in)
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	180
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes