

Suction Products, Suction Wands, Rigid Fluted Tip with Metal Tube and Plastic Handle, 21.7 Fr, SU- 12202



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

Order Code	SU-12202
------------	----------

Size	21.7 Fr
------	---------

Available Coated	no
------------------	----

Qty Per Box	10
-------------	----

Latex-Free	yes
------------	-----