

Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 21.7 Fr, SU-12402



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

Order Code	SU-12402
------------	----------

Size	21.7 Fr
------	---------

Available Coated	no
------------------	----

Effective Length (cm)	36.8
-----------------------	------

Qty Per Box	10
-------------	----

Latex-Free	yes
------------	-----