

# Suction Products, Sumps, Flexible and Malleable Sump with Wire-wound Tip, 1/4



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

Order Code	SU-20602
------------	----------

Size	1/4 in
------	--------

Available Coated	no
------------------	----

Effective Length (cm)	38
-----------------------	----

Qty Per Box	10
-------------	----

Latex-Free	yes
------------	-----