

Thoracic Catheters, Straight, 36 Fr, TC-10136



Livanova Adult Thoracic Catheters are available in straight and right angled models, and are constructed of thermo-sensitive, kink-resistant PVC material. Depth markers, in 2''' increments, aid accurate catheter placement; whereas oval drainage perforations maximize drainage performance. A blue radio-opaque stripe confirms placement upon X-ray.

ORDER CODE	TC-10136
SIZE	36 Fr
EFFECTIVE LENGTH (CM)	51,8
LATEX-FREE	yes
CONNECTOR ACCEPTANCE	3/8 in