

# Venous Return Cannulae, Dual Stage, Low Profile, Wire-reinforced Tubing, 34/46 Fr, LRD-61034



The Dual Stage Venous Return Cannula is indicated for single tube venous drainage from the right atrium and vena cava during cardiopulmonary bypass surgery. The Dual Stage Cannula has a proximal flexible polyvinyl chloride tubing segment for drainage of the right atrium. The catheter body tapers to a flexible polyvinyl chloride tubing distal segment that ends in a straight fluted distal tip for drainage of the vena cava. The proximal and distal segments of the catheter are available in wire or non-wire reinforced tubing.

<b>ORDER CODE</b>	LRD-61034
<b>SIZE</b>	34/46 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	41
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes